






# Product Catalog



# Accessories For The EGGNEST™

	Component		Description
	Name	Number	
	Mattress Extender	EN-ME-001	This is an extension from the end of the EggNest to end of the cath lab table to provide an even work surface.
	Standard Table Skirt	EN-TS-001	Customized table skirt attaches to cath lab table rail to block scatter radiation below the table.
	Table Skirt with Articulating Arm	EN-TS-002	Customized table skirt with articulating arm attaches to cath lab table rail for ideal positioning to block scatter radiation below the table.

For pricing please contact Egg Medical or your local distributor

🌐 [www.eggmedical.com](http://www.eggmedical.com) ☎ (612) 916-6616 ✉ [info@eggmedical.com](mailto:info@eggmedical.com)

To place an order please contact us at [order@eggmedical.com](mailto:order@eggmedical.com)







# Replacement Components

	Component		Description
	Name	Number	
	<b>EGGNEST™ Base Sled</b>	EN-BU-002	Carbon fiber base sled is the foundation platform to attach radiation protection components to and includes built in mattress.
	<b>Mattress Topper</b>	EN-MT-001	Independent foam mattress that's fits over and is secured to EggNest Base Sled.
	<b>Right Arm Board with Integrated Flip Shield</b>	EN-AB-R01	Carbon fiber arm board attaches EggNest platform to provide support for radial procedures and incorporates flip shields to block scatter radiation.
	<b>Left Arm Board with Integrated Flip Shield</b>	EN-AB-L01	Carbon fiber arm board attaches EggNest platform to provide support for radial procedures and incorporates flip shields to block scatter radiation.
	<b>Arm Gel Pad</b>	EN-GP-001	Gel Pad rests on the arm board to protect the patient's arm and elbow from pressure points.
	<b>Right Side Hanging Flex Shield</b>	EN-SS-R01	Hanging shielding system on right side that is supported by internal stays to allow movement of the C-arm. It attaches to the EggNest Base Sled rail.

For pricing please contact **Egg Medical** or your local distributor

🌐 [www.eggmedical.com](http://www.eggmedical.com) ☎ (612) 916-6616 ✉ [info@eggmedical.com](mailto:info@eggmedical.com)




To place an order please contact us at [order@eggmedical.com](mailto:order@eggmedical.com)

	Component		Description
	Name	Number	
	Left Side Hanging Flex Shield	EN-SS-L01	Hanging shielding system on right side that is supported by internal stays to allow movement of the C-arm. It attaches to the EggNest Base Sled rail.
	Right Head Hanging Flex Shield	EN-HS-R01	Hanging Shielding system below the table at the right side head and corner. It attaches to the rail that wraps around the head of the EggNest base platform.
	Left Head Hanging Flex Shield	EN-HS-L01	Hanging Shielding system below the table at the right side head and corner. It attaches to the rail that wraps around the head of the EggNest base platform.
	Standard Workbench with Perispherical Shield	EN-WB-001	Attachment that serves as a workbench for radial procedures and incorporates flip shielding to block scatter radiation from the lower extremity.
	Pillow Flip Shield	EN-PS-001	Shield at the head of the table that attaches to the rail and will flip up behind the pillow to block scatter radiation.
	Clear Spot Shield	EN-CS-001	Transparent Acrylic Shield on flexible gooseneck post. Positioned at the groin this device provides additional pelvic scatter radiation protection.

For pricing please contact Egg Medical or your local distributor

🌐 [www.eggmedical.com](http://www.eggmedical.com) ☎ (612) 916-6616 ✉ [info@eggmedical.com](mailto:info@eggmedical.com)

To place an order please contact us at [order@eggmedical.com](mailto:order@eggmedical.com)

	Component		Description
	Name	Number	
	Clear Spot Shield Sterile Drape	EN-SD-001	Transparent Sterile Drape that is designed to fit easily over and securely to the Clear Spot Shield on the EggNest System.
	Left Head Corner Flip Shield	EN-CS-L01	Carbon fiber shield which fits over top mattress to provide increased protection from scatter radiation at head of table.
	Right Head Corner Flip Shield	EN-CS-R01	Carbon fiber shield which fits over top mattress to provide increased protection from scatter radiation at head of table.

For pricing please contact Egg Medical or your local distributor

🌐 [www.eggmedical.com](http://www.eggmedical.com) ☎ (612) 916-6616 ✉ [info@eggmedical.com](mailto:info@eggmedical.com)

To place an order please contact us at [order@eggmedical.com](mailto:order@eggmedical.com)