



Product Catalog



Replacement Components

Name	Number	Weight		Dimensions		Description
		lb	kg	in	cm	
BASE SLED	EN-BU-002	38	17.2	83 x 29 x 2.5	210.8 x 73.7 x 6.4	Carbon-fiber Base Sled serves as a sturdy foundation platform where radiation protection components attach. This sled comes equipped with a built-in CPR board.
MATTRESS TOPPER	EN-MT-001	8.5	3.9	83 x 19.5 x 2	210.8 x 49.5 x 3.8	For maximum patient comfort, the 2" foam mattress topper is secured on top of the Base Sled.
RIGHT JAM BOARD UPPER FLIP SHIELD	EN-JU-R01	3	1.4	14 x 8.5 x 12	35.6 x 21.6 x 30.5	Carbon-fiber jam board flip shield slides under and fits against the mattress topper on the upper right side to provide protection above the table (removable for radial procedures).
LEFT JAM BOARD UPPER FLIP SHIELD	EN-JU-L01	6	2.7	26 x 10.5 x 12	66 x 26.7 x 30.5	Carbon-fiber jam board flip shields slides under and fits against the mattress topper on the upper left side to provide protection above the table (removable for radial procedures).

Replacement Components


Name	Number	Weight		Dimensions		Description
		lb	kg	in	cm	
ARM GEL PAD	EN-GP-L01	1.5	0.7	23.5 x 5 x 1	59.7 x 12.7 x 2.5	Gel Pad rests on the Arm Board to protect the patient's arm and elbow from pressure points.
RIGHT SIDE HANGING RETRACTABLE FLEX SHIELD	EN-RS-R01	9.5	4.3	40 x 15	101.6 x 38.1	Hanging Retractable Shielding System on the right side is supported by internal stays, with an adjustable skirt top allowing for greater C-arm range of motion when steep images are needed (collapses to 12"). It seamlessly attaches to the Base Sled rail and protects from scatter radiation below the table on right side.
LEFT SIDE HANGING RETRACTABLE FLEX SHIELD	EN-RS-L01	9.5	4.3	40 x 15	101.6 x 38.1	Hanging Retractable Shielding System on the left side is supported by internal stays, with an adjustable skirt top allowing for greater C-arm range of motion when steep images are needed (collapses to 12"). It seamlessly attaches to the Base Sled rail and protects from scatter radiation below the table on left side.
HEAD HANGING FLEX SHIELD	EN-HS-001	8.5	3.9	15 x 66	38.1 x 167.6	Hanging Shielding System protects from scatter radiation below the table at the head. It is securely attached to the rail that wraps around the Base Sled head.

Replacement Components

Name	Number	Weight		Dimensions		Description
		lb	kg	in	cm	
HEAD BOX	EN-HB-001	10	4.5	20.5 x 14.5 x 10	52.1 x 36.8 x 25.4	Carbon-fiber Head Box fits on head of mattress to provide increased protection from scatter radiation above the table.
STANDARD WORKBENCH	varies upon purchase					Carbon-fiber Flip Shield doubles as shielding and a workbench. Incorporates flip shielding to block scatter radiation from the lower extremity. Manufactured in various configurations to accommodate different physician/technician workflows. All configurations weight about 6.5 lbs.
RIGHT SLIDING JAM BOARD HIP FLIP SHIELD	EN-JH-R01	8.5	3.9	48 x 6.5	121.9 x 16.5	Carbon-fiber Sliding Flip Shield slides under and fits against the mattress topper on the right side to provide radiation protection above the table at the hip of the patient. Sliding allows for adjustments, can compress to 19.5".
LEFT SLIDING JAM BOARD HIP FLIP SHIELD	EN-JH-L01					Carbon-fiber Sliding Flip Shield slides under and fits against the mattress topper on the left side to provide radiation protection above the table at the hip of the patient. Sliding allows for adjustments, can compress to 19.5".

Replacement Components

Name	Number	Weight		Dimensions		Description
		lb	kg	in	cm	
CLEAR SPOT SHIELD	EN-CS-001	12	5.4	17 x 14 x 0.5	43.18 x 35.6	Transparent Acrylic Shield with adjustable arm post provides manual positional protection from the C-arm above the table.
CLEAR SPOT SHIELD STERILE DRAPE	EN-SD-001	0		Covers Clear Spot Shield		Transparent Sterile Drape that is designed to fit securely and easily over the Clear Spot Shield.



Name	Number	Weight		Dimensions		Description
		lb	kg	in	cm	
MATTRESS EXTENDER	EN-ME-001	R: 4 L: 5.5	R: 1.8 L: 2.5	20 x 24 x 4 20 x 36 x 4	50.8 x 61 x 10.2 50.8 x 91.4 x 10.2	This mattress is designed to extend from the end of the Base Sled and Mattress Topper to the end of the cath lab table, creating a seamless and even work surface. Comes in regular or long sizes.
EGG PAD	EN-EP-001	3	1.4	21.5 x 15.5	54.61 x 39.4	Reusable radiation pad placed on patient to protect from patient scatter. Effective replacement of RadPad.
EGG PAD STERILE DRAPE	EN-PD-001	0		Covers Egg Pad		Sterile Drape that is designed to fit securely and easily over the Egg Pad.
FLEXIBLE ANESTHESIA HOOP	EN-AH-001	2.25	1.0	N/A		This anesthesia hoop shapes and molds easily over the patient and table to allow for optimal workflow and drape positioning.



ORDER YOUR COMPONENTS TODAY!

 eggmedical.com

 info@eggmedical.com

 (612) 916-6616

