



Product Catalog



Replacement Components

Name	Number	Weight		Dimensions		Description
		lb	kg	in	cm	
BASE SLED	EN-BU-002	38	17.2	83 x 29 x 2.5	210.8 x 73.7 x 6.4	Carbon-fiber Base Sled serves as a sturdy foundation platform where radiation protection components attach. This sled comes equipped with a built-in CPR board.
MATTRESS TOPPER	EN-MT-001	8.5	3.9	83 x 19.5 x 1.5	210.8 x 49.5 x 3.8	For maximum patient comfort, the 1.5" foam mattress topper is secured on top of the Base Sled.
RIGHT ARM BOARD WITH INTEGRATED FLIP SHIELD	EN-AB-R01	4.5	2.0	38 x 6.5	96.5 x 16.5	Carbon-fiber Arm Board attaches to the right side of the Base Sled to offer reliable support for radial procedures. The Arm Board to equipped with innovated flip shields, effectively blocking scatter radiation and ensuring safety.
LEFT ARM BOARD WITH INTEGRATED FLIP SHIELD	EN-AB-L01					Carbon-fiber Arm Board attaches to the left side of the Base Sled to offer reliable support for radial procedures. The Arm Board to equipped with innovated flip shields, effectively blocking scatter radiation and ensuring safety.

Replacement Components

Name	Number	Weight		Dimensions		Description
		lb	kg	in	cm	
ARM GEL PAD	EN-GP-L01	1.5	0.7	23.5 x 5 x 1	59.7 x 12.7 x 2.5	Gel Pad rests on the Arm Board to protect the patient's arm and elbow from pressure points.
RIGHT SIDE HANGING FLEX SHIELD	EN-SS-R01	14	6.4	25.5 x 47	64.8 x 119.4	<p>Hanging Shielding System on the right side is supported by internal stays, allowing smooth movement of the C-arm without compromising safety. It seamlessly attaches to the Base Sled rail and protects from scatter radiation below the table on right side.</p>
LEFT SIDE HANGING FLEX SHIELD	EN-SS-L01					<p>Hanging Shielding System on the left side is supported by internal stays, allowing smooth movement of the C-arm without compromising safety. It seamlessly attaches to the Base Sled rail and protects from scatter radiation below the table on left side.</p>
HEAD HANGING FLEX SHIELD	EN-HS-001	8.5	3.9	15 x 66	38.1 x 167.6	Hanging Shielding System protects from scatter radiation below the table at the head. It is securely attached to the rail that wraps around the Base Sled head.

Replacement Components

Name	Number	Weight		Dimensions		Description	
		lb	kg	in	cm		
HEAD BOX	EN-HB-001	10	4.5	20.5 x 14.5 x 10	52.1 x 36.8 x 25.4	Carbon-fiber Head Box shield fits on top of the mattress to provide increased protection from scatter radiation above the table at the head.	
STANDARD WORKBENCH	varies upon purchase				varies upon purchase		Carbon-fiber Flip Shield doubles as shielding and a workbench. Incorporates flip shielding to block scatter radiation from the lower extremity. Manufactured in various configurations to accommodate different physician-technician workflows. All configurations typically 6.5 lbs.
RIGHT HIP FLIP SHIELD	EN-HF-R01	2	0.9	8 x 20	20.3 x 50.8	Carbon-fiber Flip Shield attaches to the right side of the Base Sled rail to provide radiation protection above the table at the hip of the patient.	
LEFT HIP FLIP SHIELD	EN-HF-L01					Carbon-fiber Flip Shield attaches to the left side of the Base Sled rail to provide radiation protection above the table at the hip of the patient.	

Replacement Components

Name	Number	Weight		Dimensions		Description
		lb	kg	in	cm	
CLEAR SPOT SHIELD	EN-CS-001	12	5.4	17 x 14 x 0.5	43.18 x 35.6	Transparent Acrylic Shield with adjustable arm post provides manual positional protection from the C-arm above the table.
CLEAR SPOT SHIELD STERILE DRAPE	EN-SD-001	0		Covers Clear Spot Shield		Transparent Sterile Drape that is designed to fit securely and easily over the Clear Spot Shield.

Name	Number	Weight		Dimensions		Description
		lb	kg	in	cm	
MATTRESS EXTENDER	EN-ME-001	R: 4 L: 5.5	R: 1.8 L: 2.5	20 x 24 x 4 20 x 36 x 4	50.8 x 61 x 10.2 50.8 x 91.4 x 10.2	This mattress is designed to extend from the end of the Base Sled and Mattress Topper to the end of the operating table, creating a seamless and even work surface. Comes in regular or long sizes.
EGG PAD	EN-EP-001	3	1.4	21.5 x 15.5	54.61 x 39.4	Reusable radiation pad placed on patient to protect from patient scatter. Effective replacement of RadPad.
EGG PAD STERILE DRAPE	EN-PD-001	0		Covers Egg Pad		Sterile Drape that is designed to fit securely and easily over the Egg Pad.
FLEXIBLE ANESTHESIA HOOP	EN-AH-001	2.25	1.0	N/A		This anesthesia hoop shapes and molds easily over the patient and table to allow for optimal workflow and drape positioning.




ORDER YOUR COMPONENTS TODAY!

 eggmedical.com

 info@eggmedical.com

 (612) 916-6616



EGG
MEDICAL